

# Contents

<b>1</b>	<b>User Manual</b>	<b>2</b>
1.1	Modelling MR signals in the presence of nervous tissue injuries	2
1.2	Detection of nervous tissue injuries . . . . .	5

# 1 User Manual

**Technical requirements:** Python 3.5, Tensorflow 1.2

**General instructions:** The programs are contained in python scripts and are intended to be executed from the terminal. To execute a program, navigate to the folder that contains it or add the directory to your PATH and follow the instructions from the appropriate section below.

## 1.1 Modelling MR signals in the presence of nervous tissue injuries

**Running the simulation** The python script that contains the simulation program is `ANNdatasynt_run.py`. To run the simulation, use the command:

---

```
python ANNdatasynt_run.py runid neuron_no schemepath -vis  
bfloatpath -syn plypath
```

---

*runid* This parameter must be specified, and should be the filename (without extension) of the Tensorflow checkpoint that contains the parameters for the trained ANN.

*neuron\_no* This parameter must be specified, and should match the number of neurons that the ANN in the checkpoint uses.

*schemepath* This parameter must be specified, and should contain the path to the directory where the experimental setup (scheme file) is located.

*plypath* This parameter should be specified if a straight synthesis should be run. `plypath` should contain the directory where the .ply files of the geometry

models are located. Either this parameter or `bfloatpath` must be specified, but not both simultaneously.

*bfloatpath* This parameter should be specified, if the synthesis should be compared to another set of data. `bfloatpath` should contain the directory where the data set that should be the basis of comparison are located.

**Example:** A sample command that runs the simulation as described in this paper is:

---

```
python ANNdatasynt_run.py synth500dir30str_90neurons 90
    simulations/20170730_grad.scheme -syn
    simulations/1.5r_all/ply/
```

---

**Output:** The synthesized diffusion signals are saved as '`nn_datasynt...`'.txt in the folder of the python script.

### Training a new ANN DWI model

The ANN can also be retrained to learn the prediction based on another experimental setup. The python script that contains the detection program is `ANNdetection_train.py`. To train a new network, use the command:

---

```
python ANNdetection\_train.py schemepath bfloatpath neuron_no -i
    iterations -l learningrate -s scale -pf
```

---

*schemepath* This parameter must be specified, and should contain the path to the directory where the experimental setup (scheme file) is located.

*bfloatpath* This parameter must be specified, and should contain the path to the directory where the diffusion signals are located. The file should be in a floating point binary format (.Bfloat).

*neuron\_no* This parameter must be specified, and should specify the number of neurons the ANN should use in its hidden layer

*iterations* The default number of iterations is 100. To specify a different number, an integer can be specified here preceded by *-i*. Note that the actual number of iterations may be less if early stoppage is prompted.

*learningrate* The default number of iterations is 0.001. To specify a different number, a float can be specified here preceded by *-l*.

*scale* With this parameter, the ANN can be trained on a random subset of the provided data set. The specified value should be between 0 and 1 and the default scale is 1.0, which takes the full data set.

*-pf* During training, the prediction is plotted every *iterations/20* steps. Specifying this parameter will turn off plotting. It is recommended to turn off plotting if training should run in the background.

**Example:** A sample command that runs the training with the reduced set as described in this paper is:

---

```
python ANNdatasynth_train.py simulations/20170730_grad.scheme  
simulations/1.5r_all/ 90
```

---

**Output:** The checkpoints of the trained network are saved in the `checkpoint` folder.

## 1.2 Detection of nervous tissue injuries

### Running the detection

The python script that contains the detection program is `ANNdetection_run.py`. To run the detection, use the command:

---

```
python ANNdetection_run.py runid neuron_no schemepath bfloatpath
--snr snr_value
```

---

*runid* This parameter must be specified, and should be the filename (without extension) of the Tensorflow checkpoint that contains the parameters for the trained ANN.

*neuron\_no* This parameter must be specified, and should match the number of neurons that the ANN in the checkpoint uses.

*schemepath* This parameter must be specified, and should contain the path to the directory where the experimental setup (scheme file) is located.

*bfloatpath* This parameter must be specified, and should contain the path to the directory where the diffusion signals are located. The file should be in a floating point binary format (.Bfloat).

*snr\_value* This parameter is optional, and can contain a signal-to-noise ratio (float). Specifying this parameter will add noise to the diffusion signals in *bfloatpath*. When specified, the value needs to be preceded by `--snr`.

**Example:** A sample command that runs the detection with the reduced set as described in this paper with a SNR of 10 is:

---

```
python ANNdetection_run.py detect10dir5str_10neurons 10
simulations/10dir5str.scheme simulations/10dir5str/ --snr 10
```

---

**Output:** The output of the detection is shown on screen.

### Training a detection ANN

The ANN can also be retrained to learn the detection based on another experimental setup. The python script that contains the detection program is `ANNdetection_run.py`. To train a new network, use the command:

---

```
python ANNdetection\_train.py schemepath bfloatpath neuron_no -i
iterations -l learningrate -nostp
```

---

*schemepath* This parameter must be specified, and should contain the path to the directory where the experimental setup (scheme file) is located.

*bfloatpath* This parameter must be specified, and should contain the path to the directory where the diffusion signals are located. The file should be in a floating point binary format (.Bfloat).

*neuron\_no* This parameter must be specified, and should specify the number of neurons the ANN should use in its hidden layer

*iterations* The default number of iterations is 100. To specify a different number, an integer can be specified here preceded by `-i`. Note that the actual number of iterations may be less if early stoppage is prompted.

*learningrate* The default number of iterations is 0.001. To specify a different number, a float can be specified here preceded by -1.

*-nostp* The default setting for early stoppage is true. Early stoppage can be turned off by adding this optional parameter.

**Example:** A sample command that runs the training with the reduced set as described in this paper is:

---

```
python ANNdetection_train.py simulations/10dir5str.scheme  
simulations/10dir5str/ 10 -i 100 -l 0.01
```

---

**Output:** The checkpoints of the trained network are saved in the `checkpoints_class` folder.